


# Chapter 4.1

# NUMERICAL

Louis J. Farrugia  
Dept. of Chemistry  
University of Glasgow  
G4 9DS, Scotland, U.K.  
email: [louis@chem.gla.ac.uk](mailto:louis@chem.gla.ac.uk)

These methods provide mathematically exact absorption corrections, providing of course that the input data are correct. The first two methods require that the crystal faces have been indexed, and the distances of the crystal faces from the centre of rotation have been measured. This can be very difficult to do accurately.


The reflection file that is actually used for input (though indicated as *name.HKL* in the following descriptions) is the file selected through the "Set HKL File" menu item. If more than one SHELX style reflection file is found, the following window is displayed :



This window indicates whether the reflection file contains direction cosines. For many absorption corrections, these are essential. The Serial data processing programs XCAD4/PROFIT will normally include direction cosines in their output reflection files, but this will not normally be the case for Area detector processed data.

The auxiliary program XtalView is provided as a useful check that the proposed set of indices and distances accurately describe the crystal. XtalView reads the file FACES.DEF (section 4.1.6) and allows the user to edit it interactively. This file is read by the programs ANALYTICAL and GAUSSIAN.

## 4.1.1 - ANALYTICAL Absorption Correction

**Synopsis**

The procedure of de Meulenaar, Tompa and Alcock as implemented in PLATON is used. The required data are the indices of the faces (which must make up a closed polyhedron) and the distance of the centre of rotation of the polyhedron to each face (in mm). The required files are name.INS and name.HKL. The latter file **MUST** contain direction cosines corresponding to the cell in the name.INS file. The direction cosines are removed in the output file.

**Files read**

name.INS (for unit cell information only).

name.HKL - uncorrected reflection file.

FACES.DEF - face indices and distances. This file may also be created interactively by the program.


**Files written**

ANALYTICAL.HKL - absorption corrected reflection file

ABSORB.LST - PLATON list file giving a list of the applied corrections

ABSORB.CIF - CIF containing relevant information

## 4.1.2 - GAUSSIAN Absorption Correction

**Synopsis**

The method of Gaussian integration as implemented in PLATON is used. The required data are the indices of the faces (which must make up a closed polyhedron) and the distance from the centre of rotation of the polyhedron to each face (in mm). The required files are name.INS and name.HKL. The latter file **MUST** contain direction cosines corresponding to the cell in the name.INS file. The direction cosines are removed in the output file.

**Files read**

name.INS (for unit cell information only)

name.HKL - uncorrected reflection file

FACES.DEF - face indices and distances. This file may also be created interactively by the program.


**Files written**

GAUSSIAN.HKL - absorption corrected reflection file

ABSORB.LST - PLATON list file giving a list of the applied corrections

ABSORB.CIF - CIF containing relevant information

## 4.1.3 - SPHERICAL Absorption Correction

**Synopsis**

This provides a correction for spherical crystals, effectively the theta dependence of the absorption surface. If  $\mu R$  exceeds 10.0 the program is aborted, as the errors will be too large. The only required input is the mean spherical radius of the crystal in mm. The correction is calculated by interpolation of data from International Tables. If direction cosines are present in the input reflection file, these will be echoed to the output file.


**Files read**

name.INS (for unit cell information only)  
name.HKL - uncorrected reflection file

**Files written**

SPHERICAL.HKL - absorption corrected reflection file  
ABSORB.CIF - CIF containing relevant information

## 4.1.4 - CYLINDRICAL Absorption Correction

**Synopsis**

This provides a correction suitable for cylindrical or needle crystals, where the crystal is mounted along the cylinder axis. If  $\mu R$  exceeds 10.0 the program is aborted, as the errors will be too large. The only required input is the mean cylindrical radius of the crystal in mm. The correction is calculated by interpolation of data from International Tables. If direction cosines are present in the input reflection file, these will be echoed to the output file.


**Files read**

name.INS (for unit cell information only)  
name.HKL - uncorrected reflection file

**Files written**

CYLINDRICAL.HKL - absorption corrected reflection file  
ABSORB.CIF - CIF containing relevant information

## 4.1.5 - NEEDLE Absorption Correction

**Synopsis**

This provides a CAD4 specific correction suitable for cylindrical or needle crystals, where the crystal is mounted along the spindle axis. It will be useful for crystals where the needle axis is very much longer than the X-ray beam width. A GUI is presented for further information and user input. If direction cosines are present in the input reflection file, these will be echoed to the output file.

**Files read**

name.INS (for unit cell information only)  
name.HKL - uncorrected reflection file

**Files written**

NEEDLE.HKL - absorption corrected reflection file  
ABSORB.CIF - CIF containing relevant information

#### 4.1.6 The FACES.DEF File

This file contains the face index and distance information about the crystal. It is a SHELX style free format text file with four-character data identifiers. It may either be created manually with a text-editor, created by the ANALYTICAL or GAUSSIAN absorption correction programs, or by the visualisation program XtalView. The format is as follows


```
TITL      test crystal shape for NaWO4.2H2O
CELL      8.4797 10.5930 13.8527  90.0  90.0  90.0
FACE      0  0  1  0.0387
FACE      0  0 -1  0.0387
FACE      1  0  0  0.1225
FACE     -1  0  0  0.1225
FACE      0  1  0  0.1375
FACE      0 -1  0  0.1375
FACE      0  1  2  0.0825
FACE      0 -1 -2  0.0825
FACE      0  1 -2  0.0825
FACE      0 -1  2  0.0825
FACE      1  0  1  0.1050
FACE     -1  0 -1  0.1050
FACE     -1  0  1  0.1050
FACE      1  0 -1  0.1050
ABSC      18.669
GRID      5
```

The TITL card is optional if you want to put a title in the picture

The GRID entry allows the user to modify the default setting of the %age accuracy required (5%). It is only used for the GAUSSIAN correction in the Area Detector programs



## 4.1.7 The XTALVIEW program



This program allows the user to visualise the crystal faces in the file FACES.DEF and also to edit the dimensions of add new faces. The viewpoint may be changed interactively by mouse dragging and the reciprocal space axial directions are also rotated